

**NITREX**

**PRODUCT INFORMATION**

<b>Item Number</b>	<b>EGN08</b>
<b>Product Colour</b>	<b>Blue</b>
<b>Size Range</b>	<b>S-XL</b>
<b>Fingertip Thickness</b>	<b>0.15 mm</b>
<b>Palm Thickness</b>	<b>0.10 mm</b>
<b>Typical Weight</b>	<b>6.50 g</b>
<b>Shelf Life</b>	<b>3 Years</b>

**FEATURES**

**Cuff:** Long cuff - 30cm glove for extra protection  
Beaded cuff to aid donning

**Grip:** Textured finish to the fingers to enhance grip and handling

**Sensitivity:** Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Powder Content:** Powder Free, reducing the risk of dust contamination

**Suitability:** Suitable for use with cytotoxic/chemotherapy drugs

**Silicone Free:** Silicone residues from gloves can interfere with painting, coating, or adhesive bonding processes, potentially causing defects. These silicone-free gloves help minimise that risk, ensuring clean application surfaces.

**Food Safe:** Suitable for food contact; all food types. Please refer to our online Food Compatibility Chart for full details

**SECTORS**

Automotive • Construction and Utility • Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Medical, Vet, and Dental



**STANDARDS**

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

**NITREX**

**CAUTION**

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

**PACKING**

100 gloves per dispenser  
10 dispensers per case  
40 cases per pallet

**INNER PACKAGING**



**STORAGE INSTRUCTIONS**

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

**DISPOSAL INSTRUCTIONS**

Dispose of the product responsibly in accordance with local regulations.

Size	Inner Packaging Barcode
Small	5060040793169
Medium	5060040793176
Large	5060040793183
XL	5060040793190

Inner Packaging	
Inner dimensions (LxWxH)	260 x 125 x 75mm

Size	Outer Packaging Barcode
Small	5060040793213
Medium	5060040793220
Large	5060040793237
XL	5060040793244

Outer Packaging	
Outer Dimensions (LxWxH)	385 x 262 x 272mm

**PLEASE NOTE:**

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555